


FOR U.S. INSTALLATION ONLY.

FOR ANCHORAGE REFER TO  
DSHPD OPA-2919-10




Weight 80 kg. - 176 lbs.

Sanitary installation		MEIKO	In the building
B	Cold Water	1/2" NPT "Y" Strainer	Inlet pipe 3/4". Stop-cocks, 41°-77°F Flow: flow pressure >4,76 Gal/Minute., >14,5 PSI
C	Hot water	1/2" NPT "Y" Strainer	Inlet pipe 3/4". Stop-cocks, 110°-140°F Flow: flow pressure >4,76 Gal/Minute., >14,5 PSI
D	Low pressure steam	R 1/2" NPT, Stop-cocks, "Y" Strainer	Ød: max. 7.3 PSI with steam trap and return Steam lines must be cleaned of debris
E	Waste Outlet TopLine20	3" "P" Trap	Waste Pipe 3" DFU 4
F	Vent TopLine20	2" Vent Connection	Local Area Vent 2" 12.4 CFM

Electrical installation		
A	1 PH ~ 208/230V, 60Hz, 0.8KW	Wall socket and main fuse. Protection: 2x15A. Protective measures as well as connection of the equi-potential bonding must be carried out according to IEC 364-5-54 and conform to the local and national power utilities regulations.

Connection dimensions are subject to tolerances.  
It is neither allowed to give this drawing without our permissions to a third party for information or reproduction nor may it be used for competition purposes. All rights are reserved! Subject to technical changes.

 <b>MEIKO USA INC.</b> 1349 Heil Quaker Blvd. LaVergne TN 37086 UNITED STATES OF AMERICA TEL. (615)399-6600; FAX: (615)399-6600			
	Alter.	Drawn	Checked

Reference	Type <b>TOPLINE20 DT-NDD</b>		
	Thermal disinfection, External Steam Supply		
	Proj. no.		
	Drawing no.		Index

Scale: 1:20	Drawn: db	Date: May.10.2011	Checked: db	Date: May.10.2011
-------------	-----------	-------------------	-------------	-------------------