



Sanitary Installation	MEIKO	In the building
B Cold Water	1/2" NPT Dirt collector	Inlet pipe 3/4", Stop-cocks, 41°-77°F Flow: Flow pressure 14,76 Gal/Minute, 14,5 PSI
C Hot water	1/2" NPT Dirt collector	Inlet pipe 3/4", Stop-cocks, 113°-140°F Flow: Flow pressure 14,76 Gal/Minute, 14,5 PSI
E Drain	Drain trap TopLine40 Dia 4 11/32"	Wall sleeve 4"

Electrical Installation

A 3 PE ~ 208/230V, 60Hz, 4,2/4,6kW

Protection class: Equipment Class 1, Pollution Degree 2, Installation Category II, Main switch,

Wall socket and main fuse. Protection: 3x15A. Protective measures as well as connection of the equal-potential bonding must be carried out according to IEC 364-5-54 and conform to the local and national power utilities regulations.

Connection dimensions are subject to tolerances. It is neither allowed to give this drawing without our permissions to a third party for information or reproduction, nor may it be used for competition purposes. All rights are reserved! Subject to technical changes.

MEIKO USA INC.
 1349 Heil Quaker Blvd, LaVerne TN 37086
 UNITED STATES OF AMERICA
 TEL: (615)399-6600 FAX: (615)399-6600

C	rm	Apr.08.2009	rm	Apr.08.2009
B	rm	Feb.19.2008	rm	Feb.19.2008
A	rm	Nov.20.2007	rm	Nov.20.2007

Alter: Drawn Checked

Reference

Type: TOPLINE40 DT-ND/V
 Thermal disinfection, decentralized

Scale:	1:20	Drawn:	rm	Date:	Jul.31.2007	Checked:	rm	Date:	Jul.31.2007
Proj. no.:	Drawing no. MTP40US-S328454				C		Index		